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

*In Vivo Assessment of Bacteria Infection Clearance of an Acellular Dermal Matrix and a Synthetic Mesh. Rosines , E; Qin, X; Chen, S. A. Symposium on Advanced Wound Care and Wound Healing Society, 27th Annual Mtg. Abr 23-28, 2014.

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
¹Barber and Aziz-Jacobo, *Arthroscopy*, 2009; 25: 1233-1239. Selected data was derived from Figure 3 (Suture Pull-out) and Table 1 (Ultimate Tensile Strength) of this reference. The two studies were performed at different points in time; however, the exact same methods, fixtures, material testing machine, and facility were used for both studies.

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Estudios *In Vitro/Vivo* con dermis Matrancell

- Capito, A.E., Tholpady, S.S., Agrawal, H., Drake, D.B. & Katz, A.J. "Evaluation of Host Tissue Integration, Revascularization, and Cellular Infiltration Within Various Dermal Substrates." *Annals of Plastic Surgery*. 2012 Oct; 68(5):495-500.
- Rosines E, Qin X, and Chen S. "In Vivo Assessment of Bacteria Infection Clearance of an Acellular Dermal Matrix and a Synthetic Mesh." Presentación participante en el Simposio sobre Tratamiento Avanzado de Heridas en la primavera de 2014.
- Agrawal H, Tholopady SS, Capito AE, Drake DB, Katz AJ. "Macrophage phenotypes correspond with remodeling outcomes of various acellular dermal matrices." *Journal of Regenerative Medicine*. 2012. Vol 1, No. 3, pp 51-59.
- Rosines, E. & Lin, Q. "Analysis of the Acellular Matrix, Growth Factors, and Cytokines Present in ArthroFlex®." 2012.
- Moore, M.A. "Inactivation of enveloped and non-enveloped viruses on seeded human tissues by gamma irradiation." 2012. *Cell Tissue Bank*. 13(3):401-7. doi: 10.1007/s10561-011-9266-0.
- LifeNet Health Research and Development "Analysis of DNA Residuals in ArthroFLEX Tissue" 2013 White Paper. Disponible bajo petición.
- Moore MA, Samsell B, Wallis G, Triplett, S, Chen S, Linthurst Jones A, Quiz X. "Decellularization of Human Dermis Using Non-Denaturing Anionic Detergent and Endonuclease: A Review." *Journal of Cell and Tissue Banking*.
- Alexis D. Armour, MD., MSc, Joel S. Fish, MD, MSc, Kimberly A Woodhouse, P.Eng., Ph.D., John L. Semple, MD, MSc. (Marzo, 2006). "A Comparison of Human and Porcine Acellularized Dermis: Interactions with Human Fibroblasts in Vitro". *Plastic and Reconstructive Surgery*.
- Hitesh, Agrawal, Tholpady, Sunil S., et. al. (Diciembre 2012). "Macrophage Phenotypes Correspond with Remodeling Outcomes of Various Acellular Dermal Matrices". *Open Journal of Regenerative Medicine*. (Diciembre), Vol 1. No. 3. 51-59(2012).

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- Beitzel, K., Chowaniec, D.M., McCarthy, M.B., Cote, M.P., Russell, R.P., Obopilwe, E., Imhoff, A.B., Arciero, R.A, & Mazzocca, A.D. “Stability of double-row rotator cuff repair is not adversely affected by scaffold interposition between tendon and bone.” *Am J Sports Med.* 2012 Mayo;40(5):1148-54.
- Van der Meijden, O.A., Wijdicks, C.A., Gaskill, T.R., Jansson, K.S., & Millett, P.J. “Biomechanical analysis of two tendon posterolateral rotator cuff tear repairs: extended linked repairs and augmented repairs.” *Arthroscopy.* 2013 Ene;29(1):37-45.
- Ely E, Figueroa, N, Gilot G. “Biomechanical Analysis of Gap Formation and Ultimate Tensile Failure Loads of Rotator Cuff Repairs with Extracellular Matrix Graft Augmentation.” *ORTHOPEDICS.* 2014;37(9):608-614.
- Ehsan, A., Lee, D.G., Bakker, A.J., 7 Huang, J.I. “Scapholunate ligament reconstruction using an acellular dermal matrix: a mechanical study.” *J Hand Surg Am.* 2012 Ago;37(8):1538-42. doi: 10.1016/j.jhas.2012.04.043.

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- Bullocks J. DermACELL: A novel and biocompatible acellular dermal matrix in tissue expander and implant-based breast reconstruction. *Eur J Plast Surg.* 2014;37(10):529-38.
- Jamal Bullocks, MD "The New Alternative For Breast Reconstruction Post Mastectomy; Case Study Series Part 1."
- Jamal Bullocks, MD "The New Alternative For Breast Reconstruction Post Mastectomy; Case Study Series Part 2."
- Jamal Bullocks, MD "The New Alternative For Breast Reconstruction Post Mastectomy; Case Study Series Part 3."
- Hanna K, Yu D, LeGallo R, Colen D, Drake D. Matracell processed, sterilized acellular dermal matrix in implant-based breast reconstruction: Clinical and histologic outcomes. *Surgical Forum Abstracts.* 2013;217(3S):S91.
- Ortiz JA. Poland's breast reconstruction with decellularized human dermal allograft. *Mil Med.* 2014;179(2):e249-52.
- Vashi, Christopher, (Marzo 2014). "Clinical Outcomes for Breast Cancer Patients Undergoing Mastectomy and Reconstruction with the Use of DermACELL, a Sterile, Room Temperature Acellular Dermal Matrix".
- Yu D, Hanna K, LeGallo R, Drake D. Comparison of the histological characteristics of ADM capsules to no-ADM breast capsules in ADM-assisted breast reconstruction. *Plastic Surgery Research Council.*

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- Mirzayan R, Conroy C, Sethi P. Distal Biceps Repair With Acellular Dermal Graft Augmentation. *Techniques in Shoulder & Elbow Surgery*. Sept 2015;16(3):89-92.
- Atzei A & Bertasi G. "Repair of Distal Triceps Tendon Rupture with a Human Acellular Dermal Matrix (ADM)." 2013 White Paper (68-20-097)
- Atzei A & Bertasi G. "Reconstruction of Finger Flexor Tendon and Pulley Repair with Human Acellular Dermal Matrix (ADM)." 2013 White Paper (68-20-098)
- Parks RM, Parks SM. Lateral Ankle Stabilization Using Acellular Dermal Allograft Augmentation. *J Am Podiatr Med Assoc*. 2015 Mayo; 105(3):209-217.
- Kokkalis, Z.T. Zanos, G., Weiser, R.W., & Sotereanos, D.G. "Trapezium resection with suspension and interposition arthroplasty using acellular dermal allograft for thumb carpometacarpal arthritis." *J Hand Surg Am*. 2009 Jul-Ago;34(6):1029-36.
- Lee, D.K. "A preliminary study on the effects of acellular tissue graft augmentation in acute Achilles tendon ruptures." *J FootAnkle Surg*. 2008 Ene-Feb;47(1):8-12.
- Yao CA, Ellis CV, Cohen MJ Kulber DA. Preserving the Posttrapeziectomy space with a Human Acellular Dermal Matrix Spacer: A Pilot Case Series of Patients with Thumb Carpometacarpal Joint Arthritis. *Plast Reconstr Surg Glob Open*. 2013 Oct; 1(7):e65.
- Rao BM, Kamal TT, Vafaye J, Taylor L. Surgical repair of hip abductors. A new technique using Graft jacket allograft acellular human dermal matrix.

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